


EXPERIMENT 2

Stoichiometry of a Chemical Reaction




DR. MIOY T. HUYNH

CHEMISTRY 136L


YALE UNIVERSITY

FALL 2018



CuSO₄

CHEMICAL REACTION



Q: *What is the formula of the precipitate?*

Q: *What are the values of w, x, y, and z?*

How can we determine these answers?

METHOD OUTLINE: VOLUMETRIC TITRATION

Prepare a standard solution of CuSO_4 :

*Dissolve a known accurate mass of $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$
in a known accurate volume of solution.*

+

A standardized solution of NaOH of
known concentration will be provided.

METHOD OUTLINE: VOLUMETRIC TITRATION

Take an accurate volume (a 10 mL or 5 mL graduated pipet) of **BLUE** copper sulfate solution.



Add phenolphthalein indicator. → Stir with stir bar.
Heat to ~70 °C and remove from hot plate.



While stirring: From buret, add **COLORLESS** NaOH to hot copper sulfate solution.



Stop adding NaOH when solution turns **PALE PINK**.



Moles of Cu^{2+} is known (volume \times concentration)
Moles of OH^- is also known.

→ $\text{Cu}^{2+} : \text{OH}^-$ mole ratio in precipitate can be figured out.

Use experimental $\text{Cu}^{2+} : \text{OH}^-$ mole ratio and charge balance principle to figure out formula of the precipitate.



Figure out the complete balanced reaction.
Figure out the net ionic equation.

Q: Can we test for “complete” precipitation of all Cu^{2+} from the solution?

Testing with S^{2-} :

CuS is BLACK and insoluble in water ($K_{\text{sp}} = 8 \times 10^{-37}$)

Do this testing in the fume hood!

When done, dispose right away in labeled waste container in the same fume hood.

Do NOT bring it to your work station.

H_2S is a toxic COLORLESS gas.

LAB REPORTS

Purpose Section	10 points
Results & Calculations Section (including graphs)	30 points
Answers to Questions	5 points
Datasheets	5 points

Lab report is due next week → Submit in Canvas as a PDF.

Report Pages are in Canvas.

Include pictures/scans of your datasheets to the end of your lab report.

Academic honesty

CLEAN UP YOUR WORK AREA. WASH HANDS.
CHECK WITH DR. HUYNH BEFORE YOU LEAVE.