



# EXPERIMENT 3

STANDARDIZATION OF THIOSULFATE

CHEMISTRY 134L // SPRING 2020

# *Method*

Volumetric titration requires

- An analyte
- A titrant
- A reasonably fast reaction between the two
- A satisfactory indicator or detection device

Titration can be used to determine the stoichiometry of a “new” reaction or the amount of an analyte if the stoichiometry is known.





# Method

Titration can be based on several types of reactions:

- Precipitation
- Acid-base neutralization
- Complex formation reactions
- Redox reactions

Among redox titrations, those based on  $I_2$  (or  $I_3^-$ , **triiodide**) are popular because **starch** is an excellent indicator for  $I_2/I_3^-$ .



# *What is molarity?*

In chemistry, molarity is a very useful way to express concentration.

$$\text{Molarity} = \frac{\text{Moles of solute}}{\text{Volume of solution (L)}}$$

Given any two quantities, we can calculate the third.

## *Standardize the thiosulfate solution*

Primary standards: very stable, high purity, and have known molar mass

Commercially available hydrates of  $\text{Na}_2\text{S}_2\text{O}_3$  are not primary standards because they may have *different* hydrates ( $\text{Na}_2\text{S}_2\text{O}_3 \cdot x\text{H}_2\text{O}$ ).




## *Standardize the thiosulfate solution*

Primary standards: very stable, high purity, and have known molar mass


Commercially available hydrates of  $\text{Na}_2\text{S}_2\text{O}_3$  are not primary standards because they may have *different* hydrates ( $\text{Na}_2\text{S}_2\text{O}_3 \cdot x\text{H}_2\text{O}$ ).



$\text{KIO}_3$  is a primary standard. It can be used to standardize thiosulfate solutions.



Titrate the  $\text{I}_3^-$  with the thiosulfate ( $\text{S}_2\text{O}_3^{2-}$ ) solution and **starch indicator**.



Now we can determine the concentration of thiosulfate solution used.

# Procedural Outline

Prepare a standard solution of  $\text{KIO}_3$ .

Known accurate mass in a known volume of solution.

*analytical balance*

*volumetric flask*





$$[\text{IO}_3^-] = \frac{\text{moles IO}_3^-}{\text{volume solution}} = \frac{\text{mass KIO}_3}{\text{molar mass KIO}_3} \times \frac{1}{0.1000 \text{ L}}$$

Pipet a known accurate volume of this standard  $\text{KIO}_3$  solution into an Eflask.

Also, prepare a solution of  $\text{Na}_2\text{S}_2\text{O}_3$  of a *not so accurate molarity*.





# Series of Titrations






# Series of Titrations




# Series of Titrations



## *An Example*

In an experiment, 0.3456 g  $\text{KIO}_3$  were dissolved in water and made up in a 100 mL volumetric flask. 10.00 mL of this was pipetted into a 250 mL Eflask and treated with excess KI and HCl. The resulting  $\text{I}_3^-$  was titrated with a  $\text{Na}_2\text{S}_2\text{O}_3$  solution of unknown molarity using starch as indicator. The blue to colorless end point was reached when 16.45 mL of the thiosulfate solution had been added. What is  $[\text{S}_2\text{O}_3^{2-}]$ ?



$$n_{\text{IO}_3^-} = [\text{IO}_3^-] \times V_{\text{pipetted}} = \left( \frac{0.3456 \text{ g KIO}_3}{214.001 \text{ g/mol}} \times \frac{1}{0.1000 \text{ L}} \right) \times (0.01000 \text{ L}) = 1.615_0 \times 10^{-4} \text{ mol IO}_3^-$$

$$n_{\text{I}_3^-} = 3 \times n_{\text{IO}_3^-} = 3 \times (1.615_0 \times 10^{-4} \text{ mol}) = 4.844_9 \times 10^{-4} \text{ mol I}_3^-$$

$$n_{\text{S}_2\text{O}_3^{2-}} = 2 \times n_{\text{I}_3^-} = 2 \times (4.844_9 \times 10^{-4} \text{ mol}) = 9.689_7 \times 10^{-4} \text{ mol S}_2\text{O}_3^{2-}$$

$$[\text{S}_2\text{O}_3^{2-}] = \frac{n_{\text{S}_2\text{O}_3^{2-}}}{V_{\text{buret}}} = \frac{9.689_7 \times 10^{-4} \text{ mol S}_2\text{O}_3^{2-}}{0.01645 \text{ L}} = 5.890 \times 10^{-2} \text{ M}$$

# *Notes*

1. Work independently today.
2. Remove gloves when using analytical balance.  
Ask Robert or Mioy to clean up spills inside balance.
3. Use mini-balance for sodium thiosulfate.  
Make sure it is set to grams (g)!