EXPERIMENT 6 The Triiodide Equilibrium

DR. MIOY T. HUYNH CHEMISTRY 136L YALE UNIVERSITY FALL 2018 All closed chemical reactions eventually reach their equilibrium state, where the speeds of the forward and reverse reactions are equal.

So there is no net change in concentrations at equilibrium.



For a given temperature, the numerical value of the equilibrium constant (K) for a chemical reaction is a constant independent of initial concentrations.

Equilibrium constants are dimensionless.

If the numerical value of K for a reaction is <u>large</u> Equilibrium lies far to the <u>right</u> <u>Products</u> dominate over reactants Such reactions go to *almost* 100% "completion"

The vast majority of reactions we have encountered in Chemistry 134L and 136L thus far fall in this category.

Today, we will study a reaction that does <u>not</u> fall in this category.

THE TRIIODIDE EQUILIBRIUM

 $I_2(aq) + I^-(aq) \rightleftharpoons I_3^-(aq)$

$$\mathbf{K} = \frac{[\mathbf{I}_3^-]_{eq}}{[\mathbf{I}_2]_{eq}[\mathbf{I}^-]_{eq}}$$

The objective is to determine the numerical value of K at room temperature.

What are the units of K?

To determine the numerical value of K,

we need to know the concentration of three species when equilibrium prevails.

We can titrate or use spectrophotometry, but there is a problem.

 $I_{2}(aq) + 2S_{2}O_{3}^{2-}(aq) \rightarrow 2I^{-}(aq) + S_{4}O_{6}^{2-}(aq)$ $I_{3}^{-}(aq) + 2S_{2}O_{3}^{2-}(aq) \rightarrow 3I^{-}(aq) + S_{4}O_{6}^{2-}(aq)$

If we titrate an aqueous solution containing <u>both</u> I_2 and I_3^- with thiosulfate, we will be able to calculate only the sum of $\{[I_2]_{eq,aq} + [I_3^-]_{eq,aq}\}$.

But this is not enough to figure out K!

<u>To overcome the problem</u>: We add an *immiscible* organic solvent such as <u>n-heptane</u> (and shake well).

After establishing the following two equilibria:

 $I_2(aq) + I^-(aq) \rightleftharpoons I_3^ I_2(aq) \rightleftharpoons I_2(n-heptane)$

Separate the two layers.

Titrate each layer separately (thiosulfate in buret & starch as indicator)



For the second equilibrium: $I_2(aq) \rightleftharpoons I_2(n-heptane)$

$$K' = \frac{[I_2]_{n-heptane}}{[I_2]_{aq}} = 41.8$$

If we know the numerical value of K', we can figure out $[I_2]_{eq,aq}$ $[I_2]_{eq,aq} = \frac{[I_2]_{n-heptane}}{K'}$

The numerical value of K' for the distribution of I₂ between water and several other *immiscible* liquids can be determined independently. They are available in the chemical literature.

How do we determine: $[I_3^-]_{eq,aq}$ and $[I^-]_{eq,aq}$



 $[I^-]_{eq,aq} = [I^-]_0 - [I_3^-]_{eq,aq} \quad (DUE \text{ TO STOICHIOMETRY})$

The Fe(III)-thiocyanate Equilibrium

 $Fe^{3+}(aq) + SCN^{-}(aq) \rightleftharpoons [Fe(SCN)]^{2+}(aq)$

<u>NOTES</u>

- 1. Organic layer titration: shake well by hand and vent.
- 2. Must do <u>three</u> aqueous layer titrations and <u>three</u> n-heptane layer titrations.
- 3. Accuracy of your experimental value of K will be graded.
- In part E, to measure 23 mL water and 2 mL solution, use your 10 mL measuring cylinder.

