



EXPERIMENT 5

CALORIMETRY (PART I)

DR. MIOY T. HUYNH // YALE UNIVERSITY
CHEMISTRY 134L // SPRING 2019

www.mioy.org/chem134

Calorimetry

All chemical reactions do one of two things:

- 1) Release energy to their surroundings
- 2) Absorb energy from their surroundings




More often than not, the energy involved is thermal energy (HEAT).

ΔH_{rxn} is the (thermal) energy lost/gained by a reaction under constant pressure.

EXOTHERMIC

If a reaction releases energy *into* the surroundings, the temperature of the surroundings go up.

Since the products are touching the surroundings, the temperature of the products also go up.




$$\Delta H_{\text{rxn}} = -480 \frac{\text{kJ}}{\text{mol}}$$



ENDOTHERMIC

If a reaction absorbs energy *from* the surroundings, the temperature of the surroundings go down.


Since the products are touching the surroundings, the temperature of the products also go down.




$$\Delta H_{\text{rxn}} = +49 \frac{\text{kJ}}{\text{mol}}$$

Today

Process B1



Process B2



Olden Days →

Specific heat of water was defined as

$$1 \text{ calorie/g}\cdot^{\circ}\text{C} = 4.184 \text{ J/ g}\cdot^{\circ}\text{C}$$

Nowadays →

Heat capacity of the “calorimeter + contents”
defined as the calibration factor.

$$\text{Units : J/K} = \text{J}/^{\circ}\text{C}$$

Calorimeter



Thermocouple temperature sensor (TTS)

→ read to 0.01 °C

Voltmeter

(measures voltage, V)

→ start at 0 V

Calorimeter
(resistor inside)

Stir plate
(stir bar inside)

Power supply &
Ammeter

(measures current, A)

→ start at 0 A

→ bring to ~0.5 V

Calibrating the calorimeter + contents

Run the reaction by mixing reactants, and measure the temperature change: ΔT_{rxn} .



Pass a known electrical current (Ampere, A) through a resistor (immersed in the calorimeter's contents) for a known amount of time (seconds, s) by applying a known electrical potential difference (Voltage, Volt, V).

Measure the temperature change: ΔT_{el} .



Calibration factor (CF)


$$\text{CF} = \frac{V \times A \times t}{\Delta T_{\text{el}}}$$



$$q_{\text{rxn}} = \text{CF} \times \Delta T_{\text{rxn}}$$

$$\Delta H_{\text{rxn}} = \frac{q_{\text{rxn}}}{n_{\text{LR}}}$$

Sample data



Notes

1. Change temperature precision to 0.01 °C
2. Change data collection duration to 1000 s.