




EXPERIMENT 8 (DAY 2)

Acids, Bases, & Buffers

DR. MIOY T. HUYNH
CHEMISTRY 136L
YALE UNIVERSITY
FALL 2018


ACIDS

PROTON
DONORS



BASES

PROTON
ACCEPTORS



If $K_a > 1$, the acid is **strong**
examples: HCl, HNO₃, ...

If $K_b > 1$, the base is **strong**
examples: NaOH, KOH, Ca(OH)₂ ...

If $K_a < 1$, the acid is **weak**

If $K_b < 1$, the base is **weak**

If $K_a \ll 1$, the acid is **very weak**
examples: CH₃COOH, HCOOH, ...

If $K_b \ll 1$, the base is **very weak**
examples: CH₃COO⁻, any conjugate base

MONOPROTIC AND POLYPROTIC ACIDS

Polyprotic acids have more than one ionizable proton (H^+).

Some Examples

Sulfuric acid is diprotic

Oxalic acid is diprotic

Citric acid is triprotic

Phosphoric acid is triprotic

⋮

Amino acids are polyprotic

Generally, **ACID-BASE** neutralization reactions are product favored.




In principle, when a **diprotic acid** is titrated with **NaOH** (a **monoprotic strong base**), two steep rises can be expected.



But, in practice, if the two K_a values are close to each other, then the two steep rises overlap.

The titration curve for **oxalic acid** with **NaOH** has two overlapping steps.




HOW DO WE DETERMINE EQUIVALENCE POINT #1?

In a **polyprotic acid** versus **NaOH** titration, how is the volume at the second equivalence point related to the volume at the first equivalence point?

Total volume of **NaOH** added to reach the second equivalence point
=
 $2 \times$ volume of **NaOH** added to reach the first equivalence point

Titration Curve (oxalic acid/NaOH)



BUFFERS

Consists of a weak **ACID-CONJUGATE BASE** pair.
They are essential in all biochemical processes.

BUFFERS have the ability to *resist* change in pH when a small amount of a strong **ACID** or **BASE** is added.

CALCULATING THE pH OF A BUFFER

$$[\text{H}^+]_{\text{buffer}} = K_a \times \frac{[\text{HA}]_{\text{eq}}}{[\text{A}^-]_{\text{eq}}} \approx K_a \times \frac{[\text{HA}]_{\text{initial}}}{[\text{A}^-]_{\text{initial}}}$$

THE HALF-EQUIVALENCE POINT

pK_a of a weak acid = pH at the half-equivalence point

In a weak **ACID (HA)**-strong **BASE** titration, at the half-equivalence point: one half of the initial moles of **HA** have been neutralized and converted to its **CONJUGATE BASE (A⁻)**; the other one half of **HA** moles are still there.

Thus, the solution in the beaker, at the half-equivalence point, is a **BUFFER** with equal **[HA]** and **[A⁻]**.

$$[H^+]_{\text{buffer}} \approx K_a \times \frac{[HA]_{\text{initial}}}{[A^-]_{\text{initial}}} = K_a$$

NOTES

1. The pH electrode is a delicate item, so handle it with care.
To work well, it needs to be well-immersed in the solution.
2. Plot both pH and its first derivative ($\partial\text{pH}/\partial V$).
→ Use the latter to judiciously choose NaOH increments.
3. In the vicinity of the EQUIVALENCE POINT, you must use very small increments of base (one drop at a time).
If not, you will not be able to determine the EQUIVALENCE POINT accurately.
4. No lab next week.